



Using your smartphone, scan code to see video of Splash! base replacement.

SPLASH BASE REPLACEMENT INSTRUCTIONS

For more information about the Splash!™ base replacement process and to see video demonstrating the replacement process, please visit www.splashpalliftinfo.com or scan the QR code.

REQUIRED TOOLS	CONTENTS OF SPLASH BASE REPLACEMENT PACKAGE
<p>3/4" deep socket wrench</p> <p>9/16" wrench</p> <p>Box cutter</p>	<p>Contact S.R.Smith if you are missing items.</p> <p>Hardware Packet</p> <ul style="list-style-type: none"> • Two—1/2" Nylock Nuts SS • Two—1/2" 316 SS Washer • One—3/8" 1.25 Hexhead • One—3/8" x 7/8" Flat Washer SS • One—3/8" Nylon Protective Washer • One—3/8" Lock Washer, Split S/S <p>Base Assembly</p> <p>Packing Tape</p> <p>Return Shipping Label</p> <p>Product Registration Form</p>
<p>BE SURE TO SAVE ALL PACKAGING MATERIALS; YOU WILL RE-USE THEM TO PACKAGE UP THE EXISTING BASE AND SEND BACK TO S.R.SMITH FREE OF CHARGE.</p>	

<p>1</p> <p>Carefully open package.</p> <p>Set aside the shipping label you will need for sending the base back to S.R.Smith.</p> <p>Tilt box on its side and carefully slide the replacement base out onto a non-scratch surface.</p> 	<p>2</p> <p>Using a 3/4" socket wrench, remove the nuts holding the handle in place, remove the handle, wood spacer and foam.</p> <p>Set aside.</p> 	<p>3</p> <p>Remove the battery from the LiftOperator.</p> <p>Disconnect the motor cord, actuator cord and hand control from the LiftOperator.</p> 	<p>4</p> <p>Remove the thumb nut and bolt on the back of the lift seat to remove the seat assembly.</p> 
---	--	---	--

5
This step may require two people.

Remove round cap to expose the base.



6
This step may require two people.

Have one person hold the mast steady. Using a 3/4" deep well socket wrench, remove the two nuts which connect the mast to the base.



7
This step may require two people.

Separate the mast from the base. Carefully place the mast on a non-scratch surface.



8
Attach the foam sheet and wood packaging spacer onto the base and tighten the bolt.



9
If the lift features a lift lock tab, remove the nut with a 9/16" wrench.



10
Use the handles to carefully pull the base out of the anchor.



11
Place the foam block onto the base stem and slide the assembly into the shipping carton.



12
Install the new base into the anchor. If a bolt was removed in step #9, re-attach lift lock anchor bolt.



13
This step may require two people.
Re-attach the mast assembly with the new hardware provided.

Using a 3/4" wrench, fully tighten the nuts.



14
Re-attach the battery and connect the motor cord, actuator cord and hand control.
Conduct a complete function test of the lift (refer to your pool lift owner's manual for complete details.)



15
Complete the product registration form, insert into the box with the old base.
Reseal the box with the packing tape provided.

SPLASH LIFT OWNER PRODUCT REGISTRATION
Please complete and return to us in the box with the Splash Lift you are returning to SRS Smith.

SPLASH POOL LIFT REGISTRATION

Contact Name: _____
 Contact Phone: _____
 Title: _____
 Facility Name: _____
 Street Address: _____
 City: _____ State: _____ Zip: _____
 Cell/Work Number: _____

Check back to confirm the new Splash Lift has been installed.
 Date Installed: _____

FAQ

Where is the product serial number located?
 The serial number will only be visible if the lift is located on the lower portion of the mast.

Where can I get more information?
 Please visit www.splashlift.com to find information about the assembly, replacement procedure and other important information. The replacement is free or contact our customer service team at (888) 497-9290, email customer@splashlift.com

16
For detailed shipping information, please refer to the enclosed orange paper marked 'RETURN SHIPPING INSTRUCTIONS.'



RESOURCES

Website:
www.splashpalliftinfo.com

Customer Service:
(888) 497-9290

FedEx Website:
www.fedex.com

FedEx Pick Up:
(800) 463-3339

WARRANTY INFORMATION

The Splash! base replacement assembly is covered by a two year warranty upon receipt.

We are also extending the original warranty for structural components of your lift from three years to five years from date of manufacture.

For more information, visit:

<http://srsmith.com/en-us/customer-service/warranty-information/>